# EOLife<sup>®</sup>

---

INSTRUCTION D'UTILISATION  
INSTRUCTIONS FOR USE

ARCEON  
MEDICAL

# EOLife®

---

[FR]	Instruction d'utilisation . . . . .	4
[EN]	Instructions for use . . . . .	27
[NO]	Bruksanvisning . . . . .	50
[SV]	Bruksanvisningen . . . . .	73
[FI]	Käyttöohjeet . . . . .	96
[DA]	Brugsanvisningen . . . . .	119
[DE]	Gebrauchsanweisungen . . . . .	142
[IT]	Istruzioni d'uso . . . . .	165
[ES]	Instrucciones de uso . . . . .	188
[PT]	Instruções de utilização . . . . .	211
[NL]	Gebruiksaanwijzing . . . . .	234
[KO]	사용 설명서 . . . . .	257

# TABLE DES MATIÈRES

FR

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
## TABLE OF CONTENTS

I. Symbols used on product labelling . . . . .	28
II. Acronyms . . . . .	28
III. Conditions of use . . . . .	29
1. Intended purpose . . . . .	29
2. Indication for use . . . . .	29
3. Patients. . . . .	29
3. Environment of use . . . . .	29
5. Users. . . . .	29
6. Contraindications/Restriction on use . . . . .	29
IV. Safety. . . . .	29
V. Product description . . . . .	31
VI. Equipment . . . . .	32
VII. How to use the device . . . . .	32
1. Initial product set-up. . . . .	32
2. Using EOLife® . . . . .	33
VIII. Detailed product description and operating instructions . . . . .	34
1. Using the device. . . . .	34
a) Assembly/disassembly . . . . .	34
b) EOLife® ventilation method. . . . .	35
c) Alarm conditions and alarm signals . . . . .	39
d) Standby mode . . . . .	41
2. EOLife® device set-up screen. . . . .	41
a) Set-up screen display. . . . .	41
b) Update . . . . .	42
3. Cleaning and maintenance. . . . .	42
a) Charging the battery . . . . .	42
b) Cleaning EOLife® and its accessories. . . . .	43
c) Storing EOLife®. . . . .	44
d) Update . . . . .	44
e) Disposal . . . . .	44
IX. Technical data . . . . .	45
1. Technical data relating to the device. . . . .	45
2. Technical data relating to the power supply . . . . .	46
3. Technical data relating to the FlowSense® consumable. . . . .	47
4. Technical data related to electromagnetic compatibility (EMC). . . . .	48
X. What to do in the event of an incident . . . . .	49
XI. Warranty and limits of liability . . . . .	49
XII. CE declaration of conformity for medical devices . . . . .	49

EN

## I. Symbols used on product labelling

Symbol	Definition
	Manufacturer's lot number
	Manufacturer
	Part Number
	Serial number
	Use by date
IP44	Ingress protection against: <ul style="list-style-type: none"><li>o Intrusion of solid objects larger than 1 mm</li><li>o Water splashing from any direction</li></ul>
	Do not re-use
	Do not dispose of with general waste
	CE marking
	Temperature limitation
	Humidity limitation
	Medical device
	Warning
	Do not use if package is damaged
	Quantity

## II. Acronyms

AHA	American Heart Association
BLS	Basic life support
WEEE	Waste Electrical and Electronic Equipment
ERC	European Resuscitation Council
LED	Light-Emitting Diode
CPR	Cardiopulmonary resuscitation

### III. Conditions of use

---

#### 1. Intended purpose

EOLife® is a medical device to help healthcare professionals perform pulmonary ventilation during cardiopulmonary resuscitation and adapt their technique to comply with international guidelines from the AHA and ERC.

#### 2. Indication for use

EOLife® is intended for use during cardiopulmonary resuscitation of cardiopulmonary arrest patients.

#### 3. Patients

EOLife is intended to be used on adult patients between 140 cm and 200 cm

#### 3. Environment of use

EOLife is indicated to be used in pre-hospital care, during emergency patient transport in ambulance and in hospital environment.

EOLife® is intended for use in a home healthcare environment according to EN 60601-1-2

#### 5. Users

EOLife is indicated to be used by healthcare professionals certified in Basic Life Support (BLS) or Advanced Life Support (ALS)

Users must inform himself about the safe use of the device before using it by reading the instructions for use

When the device is delivered to a department by an ARCHEON staff member or approved distributor, training should be provided in the handling and use of EOLife®

#### 6. Contraindications/Restriction on use

EOLife® is not approved for the following uses:

- In a hyperbaric chamber
- Near magnetic resonance imaging machines
- Air ambulance
- Pediatric and neonatal populations

EN

### IV. Safety

---

#### Precautions

- Before using the product, read these instructions carefully. Failure to follow the instructions can lead to improper use and expose the patient, user or persons nearby to risks.
  - Use the device as intended (see “III. Intended use”).
  - Comply with the restrictions on use (see “III. Intended use”).
  - Carefully read the safety instructions in these instructions for use.
  - Comply with all the sections in these instructions.

- The use of defective devices could compromise the ventilation of the patient and impact their condition and may also put the user at risk.
  - Only use the EOLife® device and its accessories if they do not show any visible damage.
- The EOLife® device, its battery, charger and the FlowSense® sensor must be stored and used as described in the instructions for use (see "IX. Technical data").
- We advise carrying the EOLife® device in its carrying-case or kit-bag when included in a first aid bag (see "VI. Equipment and accessories").
- The device can be safely used outdoors, as it is rainproof, splashproof and cannot be penetrated by solid objects larger than 1 mm (IP44 rating). However, it should be handled with care and kept away from powerful water jets and dust.
- Always check the battery level before using EOLife®. We advise recharging the battery after each use. Change the battery if the charge is below 20%. We also advise having a spare operational battery in case this is needed.
- Using a battery with a low charge at temperatures below 0°C greatly reduces the run time and can cause EOLife® to stop functioning after a short time. If the device is used at low temperatures, the battery should be fully charged
- Do not remove the battery when EOLife® is switched on.


## Warnings

- ⚠ Avoid using EOLife® next to or in combination with other devices as this could cause it to malfunction. If devices must be used next to each other or together, monitor both devices to check that they are functioning properly.
- ⚠ Using the device in a way other than as described by ARCHEON and using non-approved accessories could affect the patient's ventilation and health.
- ⚠ To ensure that EOLife® functions effectively you must follow the instructions provided by ARCHEON.
- ⚠ Never use the EOLife® device or its accessories in a hyperbaric chamber as this could cause an explosion. Never use the EOLife® device or its accessories in a hyperbaric chamber.
- ⚠ Any technical alterations to the device could affect the patient's ventilation and health.
  - No alterations should be made to the device or its accessories.
- ⚠ The use of accessories, transducers and cables other than those specified or supplied could increase the electromagnetic emissions or decrease the electromagnetic immunity of the device and cause it to malfunction.
- ⚠ Portable high frequency communication devices (e.g. wireless devices, antennas and their cables) in close proximity to EOLife® can interfere with the device and reduce its performance. Keep EOLife® at least 30 cm from portable high frequency communication devices.
- ⚠ Handling the battery improperly could seriously injure the patient, user and persons nearby. The battery must not be disposed of in a fire, opened, deformed or short-circuited. The battery must be protected from humidity, high temperatures and high pressure (see "IX. Technical data").
- ⚠ The FlowSense® flow sensor is for single use only. It must be used as soon as it is removed from its packaging and discarded after use. Re-using or decontaminating single-use items could contaminate the patient or have unforeseeable life-threatening consequences. It is advisable to always carry a spare FlowSense® sensor.
- ⚠ If you notice any malfunction or damage to the product or packaging, do not use EOLife® or the FlowSense® sensor and return them to the manufacturer or distributor.
- ⚠ Always clean the EOLife® electronic unit after use as instructed in "VIII. 3 b) Cleaning EOLife® and its accessories".
- ⚠ The EOLife® electronic unit, the FlowSense® sensor, the battery and the charger must be discarded at the end of their lifetime as instructed in "VIII. 3 d) Disposal".

## V. Product description

- o The device is a portable machine that includes the Eolife® electronic unit, its removable and rechargeable battery, its charger and the single-use FlowSense® sensor (the expiry date is specified on the outer packaging). Eolife® is an electro-medical device.

EN

- o Eolife® is used with the FlowSense® sensor which is positioned between a non-invasive ventilation interface (mask), invasive ventilation interface (endotracheal tube), or any type of supraglottic airway device, and a standard self-inflating bag to manually ventilate adult patients.


### Note

The ventilation interface and the self-inflating bag are accessories not supplied by ARCHEON.

### Note

The device's inlet and outlet ports comply with the international EN ISO 5356-1 standard. Eolife® is therefore compatible with any type of ventilation interface and self-inflating bag.

- o Eolife® allows the user to choose the resuscitation mode (continuous ventilation or CPR with a 30:2 compression/ventilation ratio) and the patient's height. Based on the inspiratory and expiratory flow rates measured during ventilation, Eolife® calculates the main ventilation parameters (insufflation volumes, tidal volumes, ventilation frequency and leakage) and gives real-time feedback on the ventilation performed. Eolife® also shows the target ventilation values (in accordance with the guidelines issued by the AHA and ERC international scientific committees) to guarantee adequate ventilation.

## VI. Equipment

Designation	Sales Part Number
EOLife®	A0000055
Battery	A0000051
Charger and 4 adapters	A0000029
FlowSense® (x10) (Single use component)	A0000044

The equipment can be ordered from approved distributors.

## VII. How to use the device

This section gives a short step-by-step description of how to use the EOLife® device. There is a more detailed description in the next section, "VIII. Detailed product description and operating instructions".

### 1. Initial product set-up

**Note**

The initial set-up process must be implemented only once, when you receive the EOLife® device.

**Note**

EOLife® has a touch screen that can be used while wearing gloves.

1

- Charge the battery fully  
 • See "VIII 3 a) Charging the battery"

3



- Press the physical ON/OFF button to show the summary screen

4

- Display the set-up screen


2

- Start up EOLife® by pressing the physical ON/OFF button


5

- Select the date, time and language and choose the desired interface  
 • See "VIII 2 a) Set-up screen display" and "VIII 2 b) Choice of interface"

Summary screen


Set-up screen


## 2. Using EOLife®

Assemble EOLife®, attaching the FlowSense® sensor and other devices to ventilate the patient (ventilation interface and bag)


- See “VIII 1 a) Assembly”

**1** EOLife® assembly with the FlowSense® sensor


**2** Assembly with the ventilation interface

**3** Assembly with the bag

- 4** Start up EOLife® by pressing the physical ON/OFF button

- 5** Press on the height of the adult patient in cardiopulmonary arrest
- Small: the patient's height is between 140 and 160 cm
  - Medium: the patient's height is between 160 and 180 cm
  - Tall: the patient's height is between 180 and 200 cm

Patient's height		55%
	<b>Small</b> 140 / 160 cm	
	<b>Medium</b> 160 / 180 cm	
	<b>Tall</b> 180 / 200 cm	

- 6** Choose the type of resuscitation performed

- 30:2 mode: alternating between 30 chest compressions and 2 ventilations
- Continuous mode: continuous ventilation at a frequency of 10 breaths per minute.

Ventilation Mode		55%
	<b>30 : 2</b>	
	<b>Continuous</b>	

- 7** Ventilate the patient in cardiopulmonary arrest for as long as necessary
- See “VIII 1 b) EOLife® ventilation method” and “VIII 1 c) Alarm conditions and alarm signals”

Medium						55%
<b>Vi</b>	<b>Vt</b>	<b>Freq.</b>				
420 mL	390 mL	11/min				
<b>Good Ventilation</b>						

- Information area
- Vi = insufflated volume
  - Vt = tidal volume
  - Freq. = ventilation frequency

Bar graph (real-time Vi)


Alarm signal display area

Main screen



- 8 Switch off EOLife® by pressing the ON/OFF button, then on the OFF icon or by pressing the ON/OFF button for 3 seconds.

- 9 Disconnect the ventilation bag and interface (mask or endotracheal tube) and unclip the FlowSense® sensor by pulling the tab outwards, as shown by the arrow.

EN

Switching off using the summary screen

Disconnecting EOLife®

## VIII. Detailed product description and operating instructions

EOLife® is controlled via the touch screen and the physical ON/OFF button.

### Caution

Avoid using any object that could scratch the screen.

#### 1. Using the device

##### a) Assembly/disassembly

##### Assembly

Remove the FlowSense® sensor from its packaging, using the notch to tear the antistatic bag to avoid damaging the sensor.


### Caution

Check the sensor's expiry date before using. The use by date can be found on the label affixed to the bag . Check that the bag is intact before using the FlowSense® sensor. If the FlowSense® sensor bag is torn, discard the sensor and use another one as water or contaminants may have entered the bag.

Connect the FlowSense® sensor to the EOLife® unit, taking care not to expose the printed circuit board to dust or splash it with water.


First insert the part of the FlowSense® with a notch into the connector ridge, and then press firmly on the sensor so that it comes into contact with the connector pins. You will hear a click when the connector tab holds the FlowSense® sensor in place.

Check that the FlowSense® sensor is held properly in place before starting up the EOLife® unit.

Connect the ventilation interface (mask or endotracheal tube) to the compatible end of the FlowSense® sensor. The sensor ends are shaped differently to ensure that the interface cannot be connected to the wrong end.

Connect the bag to the other end of the FlowSense® sensor. The device must be assembled as shown below:

## Disassembly

To remove the accessories from the product, take the FlowSense® sensor out of the ventilation interface (mask or endotracheal tube) by gently pulling vertically on the FlowSense® sensor. Then remove the ventilation bag at the other end in the same way.

### **Caution**

If the patient is intubated and secretions appear in the FlowSense® sensor during ventilation: remove the sensor very carefully from the endotracheal tube (or other supraglottic airway device) by holding the tube firmly with one hand and gently pulling the FlowSense® sensor upwards with the other hand to disconnect it from the tube. To prevent substances or liquids from entering the FlowSense® (patient regurgitation, patient secretion, condensation, etc.), which could affect the flow measurements and cause minor injuries to the patient, we recommend placing a filter between the ventilation interface (mask or endotracheal tube) and the FlowSense®.

## b) EOlife® ventilation method

EOlife® gives the user real-time feedback on the quality of the ventilation performed on the patient based on a cyclic analysis of inspiratory and expiratory flow rates. It must therefore be positioned proximally to be able to measure both the inspiratory and expiratory flow rates, calculate the volumes insufflated and expired, the ventilation frequency and leakage, estimate the tidal volumes, and generate user information signals and alarm signals if the ventilation parameters are unacceptable.

EOlife® also uses algorithms to filter out small airflows generated by chest compressions so as not to display anomalous volume and ventilation frequency values and to avoid measurement artefacts. During ventilation, be sure to generate enough flow to be detected by the device.


### **Warning**

When ventilating a patient with a mask, ensure it is correctly held on the patient's face during insufflation and for the entire duration of expiration so that EOlife® can measure the volume expired.

### **Note**

If used for training on a manikin, check that the manikin's airways are completely sealed so that all of the air insufflated via the mouth is fully expired. If this is not the case, EOlife® will not be able to analyse the ventilation cycles or calculate the ventilation parameters correctly.

The main screen guides the user by displaying the ventilation parameters and ventilation instructions and visual indicators showing the ideal air/oxygen volume to administer and when to administer it. The main screen also uses alarm signals to warn the user if the patient is not ventilated correctly.

1. Displays the patient height selected and quick access key to change the patient's height.
2. Displays the chosen resuscitation mode and quick access key to change the resuscitation mode.
3. Rotation icon and key to rotate the screen 180° to adjust the direction of the screen to the user's position during ventilation.  
The user can perform ventilation and change the patient's height or the resuscitation mode in this position. The screen returns to its normal position if the user:
  - Presses again on the rotation icon
  - Presses the physical ON/OFF button to show the summary screen
  - Presses the physical ON/OFF button for 3 seconds to switch the device off.
4. Displays the battery level (between 0% and 100%).  
o See "VIII. 3 a) Charging the battery".
5. On interface: displays the actual volume insufflated (Vi) for each ventilation cycle (value in mL).
6. On interface: displays the tidal volume (Vt) (value in mL):
  - o For 30:2 mode → mean tidal volume based on the two consecutive ventilation cycles.
  - o For continuous mode → actual value of tidal volume at each ventilation cycle.


**Note**  
The tidal volume is the estimated volume of air/oxygen actually reaching the patient's lungs based on both the measurement of the volume of air expired and the leakage calculated during the insufflation/expiration phases.  
The ERC and AHA recommend that an adult patient in cardiopulmonary arrest should be ventilated with a tidal volume of 6 to 8 mL/kg of theoretical ideal weight.

7. On interface: displays the trend value of the ventilation frequency (Freq.) (breaths/minute) taking into account the frequency of the last ventilation cycles performed:
  - o For the 30:2 mode → value available from the third ventilation cycle.
  - o For the continuous mode → value available from the fourth ventilation cycle.

**Note**  
The ERC and AHA recommend that an adult patient in cardiopulmonary arrest should be ventilated at a frequency of 10 breaths per minute during continuous cardiac massage.

8. Real-time indication of how to ventilate:

- o A bar graph shows whether the volume insufflated/expired is sufficient:


- When the user ventilates the patient, the bar graph shows in real time the quantity of air/oxygen that must be insufflated to the patient based on ERC and AHA guidelines.

**Note**

- The bar graph is displayed on each breath and is incremented in proportion to the volume insufflated.
  - During expiration, the bar graph is decremented in proportion to the measured volume of air/oxygen expired.
  - Take care to generate enough flow for the system to detect the insufflation and ensure that it is not mistaken for an artefact or an involuntary movement of air caused by chest compressions or handling the bag.
- o Rate display in continuous mode:

- In continuous resuscitation mode, EOLife® generates a visual countdown to show the user when to ventilate to deliver 10 breaths per minute. When the green lung is displayed, the user must ventilate the patient.

- o Display in 30:2 mode:

In 30:2 mode, the rate function is disabled as the person ventilating must wait for 30 chest compressions before delivering 2 ventilations. In this mode the ventilation frequency depends on the person (or machine) performing the chest compressions. The green lung is shown permanently to indicate that the device is waiting to detect a ventilation to display the bar graph.

9. Alarm signal display area:

- o "Good ventilation": all the ventilation parameters calculated are within the appropriate range of values as defined in the international guidelines.
- o "Insufficient volume": the volume of air insufflated to the patient is significantly lower than the values given in the international guidelines.
- o "Excessive volume": the volume of air insufflated to the patient is significantly higher than the values given in the international guidelines.
- o "Low frequency": the ventilation frequency is much significantly than the values given in the international guidelines.
- o "High frequency" (for continuous mode): the ventilation frequency is significantly higher than the values given in the international guidelines.
- o "Important leakage": the difference calculated between the volume of air/oxygen insufflated and the volume expired is too great. This may be because the mask is not positioned correctly, due to improper intubation (oesophageal intubation) or because the ventilation circuit is not properly sealed.

**Note**

See details of alarms in "VIII. 1 c) Alarm conditions and alarm signals".

## c) Alarm conditions and alarm signals

### Alarm messages related to ventilation parameters

Alarm messages may appear during ventilation when the ventilation parameters exceed the threshold values and are persistently outside the ranges given below.

The screen shows the alarms in the alarm signal display area (lower part of the screen) as a message.

The alarm messages concern:

- **The tidal volume (Vt):**

Alarm conditions	Cause		Solution
	Patient height	Tidal volume value	
Insufficient volume	Small	< 225 mL	Insufflate a volume of air within the appropriate range as defined in the international guidelines.
	Medium	< 285 mL	
	Tall	< 350 mL	
Excessive volume	Small	> 430 mL	
	Medium	> 565 mL	
	Tall	> 710 mL	

**Note**

- In continuous resuscitation mode the alarm condition is calculated based on the trend value of the volumes. If the trend value and the current value of the tidal volume are as described in the table above, the tidal volume alarm is triggered.
- In 30:2 CPR mode the alarm condition is calculated based on the mean value of the tidal volumes during two consecutive ventilation cycles.

- **The ventilatin frequency (Freq.):**

Alarm conditions	Ventilation frequency value	Solution
Low frequency (continuous mode)	No breaths for 9 seconds	Ventilate the patient at an appropriate frequency, as recommended in the international guidelines
Low frequency (30:2 mode)	No breaths for 24 seconds	
High frequency	> 15 bpm	

**Note**

- In continuous resuscitation mode the frequency is calculated based on the trend value of the ventilation frequency. If the trend value and the current value of the frequency are as described in the table above, the ventilation frequency alarm is triggered.
- In 30:2 CPR mode the alarm condition is calculated based on the actual value of the pause time.

- **Leakage:**

Alarm conditions	Leakage value	Solution
Important leakage	> 40 %	Depending on the cause of the leakage, several solutions are possible (non-exhaustive list): <ul style="list-style-type: none"> <li>• Reposition the tracheal tube if the patient is improperly intubated (selective intubation, tube not deep enough inside the trachea) or inflate the bag more if the ventilation circuit is not properly sealed</li> <li>• Reposition the mask properly on the patient's face to ensure that the ventilation circuit is properly sealed</li> <li>• Completely free the patient's airways</li> </ul>

**Note**

- In continuous resuscitation mode leakage is calculated based on the trend value of the leakages. If the trend value and the current leakage value are as described in the table above, the leakage alarm is triggered.
- In 30:2 CPR mode the alarm condition is calculated based on the mean value of the leakages during two consecutive ventilation cycles.

- **Alarm prioritisation :**

All alarm messages appear in the same display area, therefore a system is necessary to prioritise the alarms. Alarm prioritisation also enables the user to focus first on correcting the most critical problem in terms of patient safety, and avoids the additional stress of multiple simultaneous alarms.

The alarms were prioritised with clinical experts, using a structured risk management process.

Multiple alarm conditions			Priority alarm
Ventilation frequency	Tidal volume	Leakage	
Low frequency	Insufficient volume		Low frequency
Low frequency	Insufficient volume	Important leakage	Low frequency
Low frequency			Low frequency
Low frequency		Important leakage	Low frequency
Low frequency	Excessive volume		Low frequency
Low frequency	Excessive volume	Important leakage	Low frequency
	Insufficient volume		Insufficient volume
	Insufficient volume	Important leakage	Important leakage
		Important leakage	Important leakage
	Excessive volume		Excessive volume
	Excessive volume	Important leakage	Excessive volume
High frequency	Insufficient volume		High frequency
High frequency	Insufficient volume	Important leakage	High frequency
High frequency			High frequency
High frequency		Important leakage	High frequency
High frequency	Excessive volume		High frequency
High frequency	Excessive volume	Important leakage	High frequency

## Alarm message linked to the FlowSense® sensor

During ventilation, if the FlowSense® sensor is defective or not properly connected or if a foreign object on the connector interferes with the signal transmission, an alarm screen is displayed:

Alarm screen  
"FlowSense failure"

You must reconnect a functional FlowSense® sensor as described above to clear the "FlowSense failure" alarm.

### Note

The "FlowSense failure" alarm screen can only be displayed from the main screen.

## EOLife® device failure

If there is an EOLife® device failure, the red LED in the top left corner of the screen comes on. This red LED tells you when there is a device failure whether it is on or off.

### Warning

Do not use the EOLife® device if the red LED stays on.  
See "X. What to do in the event of an incident"


## d) Standby mode

If EOLife® remains in standby mode for more than 1 hour, the device will turn itself off completely.

## 2. EOLife® device set-up screen


### a) Set-up screen display

1. Start up EOLife® by pressing the physical ON/OFF button.
2. Press normally the physical ON/OFF button to show the summary screen.


Summary screen

3. Next, select the [wrench ] to display the EOlife® set-up screen.

Summary screen

Set-up screen

### b) Update

EOlife® has software that can be updated via a cable provided on request. When an update is available, you will be notified by your authorized dealer or the manufacturer.

## 3. Cleaning and maintenance


### a) Charging the battery

Recharging the EOlife® battery:




**Note**  
The battery must be removed from the EOlife® unit before recharging

1. Disconnect the battery .
2. Plug the jack on the charger cable provided into the battery.

3. Plug the charger cable to the mains socket using the appropriate adapter for your country.

Adapter name	Male adapter B (USA)	Male adapter E (Europe)	Male adapter G (UK)	Male adapter I (Australia, NZ)
Image				
Compatible sockets	Type A Type B	Type C Type E Type F	Type G	Type I

 **Warnings**

- o Only recharge the EOLife® battery using the charger provided.
- o Only use the EOLife® battery referred to on the list of accessories and equipment (part number not available).
- o If the charge is too low, do not interrupt ventilation to replace the battery. Continue ventilation without real-time feedback from EOLife®.

Always check the battery level before placing the device in the emergency bag. The battery level is shown in the top right corner of the start-up screen.

The battery icon fills as the charge increases and also changes colour:

	Green Run time ≥ 50%	Orange 50% > Run time ≥ 20%	Red 20% > Run time
Charge			

 **Note**  
Under normal conditions of use and when the operating conditions specified in "IX. Technical data" are met, the run time is around 5 hours from fully charged.

 **Caution**

When the  symbol appears on the screen, the battery is probably defective, and you should use another battery to avoid a power failure during resuscitation.

**b) Cleaning EOLife® and its accessories**

- o **FlowSense®**  
The FlowSense® sensor is a single-use component of EOLife® (see "VIII. d) Disposal").
- o **Electronic unit and removable battery:**  
The electronic unit and the removable battery are reusable EOLife® components that must be cleaned after each use to avoid cross-contamination.

This mandatory and essential step is the user's responsibility.

 **Caution**

The EOLife® device must be switched off before being cleaned.

## Cleaning

EOLife®'s electronic unit and battery are reusable. It is recommended to clean them after each use. The battery must be removed from EOLife®'s electronic unit. Before applying the cleaning solution, remove any remaining adhesive material from the device and wipe it thoroughly using a cloth moistened with water. Then use a cloth moistened with the cleaning solution and wipe all accessible surfaces of the device and its battery. No adhesive material should be visible to the naked eye under normal light on any part of the product's surface after cleaning.

If any adhesive material is still visible, repeat the cleaning process.

### Warning

Do not immerse the device. Do not pour or spill liquids onto the device. EOLife® and its battery are compatible with solutions such as "SURFANIOS Premium" and "Septalkan" containing didecylidimethylammonium chloride and isopropyl alcohol.

### Note

EOLife® has been tested for resistance to isopropyl alcohol in accordance with the requirements of EN 60601-1.

## c) Storing EOLife®

---

After cleaning the electronic unit and the removable battery, check the battery level (see "VIII. 3 a) Charging the battery").

Once charged, reconnect one of the batteries to the electronic unit and store EOLife® and its accessories in a place that meets the storage conditions (see "IX. Technical data").

## d) Update

---

EOLife® has software that can be updated via a cable provided on request. Your approved distributor or the manufacturer will notify you when an update is available.

EOLife® has a data transfer feature. To transfer the data recorded by EOLife® ask your approved distributor and/or contact the manufacturer ARCHEON.

## e) Disposal

---

At the end of the EOLife® device's lifetime (lifetime of the various parts specified in "IX. Technical data"), it should be disposed of as follows:

1. The EOLife® electronic unit, battery and charger contain electronic components and must be disposed of in an appropriate recycling facility, in accordance with European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). They must not be disposed of with household waste. Take them to an electrical and electronic equipment recycling collection point. They must be disposed of in accordance with local environmental regulations on waste disposal. For more detailed information on the handling, collection and recycling of this product, contact your environmental officer or local council.
2. Users must dispose of the FlowSense® sensor as biohazardous waste.

## IX. Technical data


### 1. Technical data relating to the device

Specifications	EOLife®
Product class according to Regulation (EU) 2017/745	Class I
Dimensions (W x H x D)	135 mm x 75 mm x 35 mm
Weight including battery	150g ± 5 g
Operating conditions	<ul style="list-style-type: none"><li>○ Temperature between 0°C and +40°C</li><li>○ Relative humidity between 15% and 95% (non-condensing)</li><li>○ Atmospheric pressure range of 620 hPa (at an altitude of 4,000 m) to 1,060 hPa (at an altitude of -500 m)</li></ul>
Transient operating conditions (maximum 20 min)	<ul style="list-style-type: none"><li>○ Temperature between -20°C and +50°C</li><li>○ Relative humidity between 15% and 90% (non-condensing)</li></ul>
Storage/transport conditions	<ul style="list-style-type: none"><li>○ Temperature between -40°C and +70°C</li><li>○ Relative humidity between 10% and 90% (non-condensing)</li></ul>
Lifetime	5 years
Power supply (nominal voltage)	3.7 V
Run time	5 hours
Classification according to EN 60601-1: <ul style="list-style-type: none"><li>○ Type of protection against electric shock</li><li>○ IP rating for protection against electric shock</li></ul>	The entire device except for the charger (electronic unit, battery and FlowSense® sensor) has been designed to meet the requirements for type BF applied parts.
IP rating for protection against solids, dust and ingress of water	IP44 (in use, i.e. when EOLife® is connected to its battery and to the FlowSense®)
Electromagnetic compatibility (EMC) according to EN 60601-1-2	The control parameters and threshold values can be obtained from the manufacturer.
Shock and vibration resistance	EN 60601-1-12 (category: resistance in an emergency vehicle)
Screen	60 x 45 mm Resolution 320 x 240 pixels

Applicable standards	EN 60601-1:2006/A1:2013/A12:2014 EN 60601-1-2:2015/A1:2021 EN 60601-1-12:2015 EN 62366-1:2015 EN 62304:2006/A1:2015 ISO 18562-1:2024 ISO 18562-2:2024 ISO 18562-3:2024
Essential performance requirements	Under normal conditions of use, EOLife® must display the bar graph in real time. The screen must not freeze.
Operating range	Vi : 50 à 2000 mL Vt : 0 à 2000 mL Freq : 0 à 60 cycles/min ± 1 cycle/min
Measurement accuracy	<p>Volume measurements are based on FlowSense® flow sensor measurements and are expressed in mL for the BTPS (body temperature and pressure, saturated). The measurement accuracies of the values displayed on the screen are as follows:</p> <ul style="list-style-type: none"> <li>• Vi (volume insufflated): ± 8.8% of the actual value measured under normal conditions of use</li> <li>• Vt (tidal volume) without leakage: ± 8.8% of the actual value measured under normal conditions of use</li> <li>• Freq (ventilation frequency): ± 1 cycle per minute</li> </ul> <p>FlowSense® data:</p> <ul style="list-style-type: none"> <li>• Flow range: ± 250 slm (standard litre per minute)</li> <li>• Dead space: &lt; 10 ml</li> <li>• Update time: 0.5 ms</li> </ul> <p>Note: Certain types of ventilation bag can affect the measurement accuracy due to their design (non-laminar outgoing air flow). A slight measurement deviation may be seen but which has no impact on compliance with regulatory requirements.</p>

## 2. Technical data relating to the power supply

Specifications	EOLife® battery
Dimensions (W x H x D)	73 mm x 40 mm x 14 mm
Weight	50 g ± 5 g
Shelf life	6 months (at a storage temperature between -10°C and +35°C) <i>Warning:</i> Battery life may be reduced if it is stored outside this temperature range. Battery life can be extended by charging the battery.
Type	Li-ion
Nominal capacity	1,900 mAh (7.03 Wh)
Nominal voltage	3.7 V
Charging voltage	12V
Charge time (0% to 95%)	Approx. 2 hours
Recommended charging temperature	+10°C to 45°C
Certifications	IEC 62133 UN38.3 UL1642

Specifications	Eolife® Charger			
Weight	206 g including accessories			
Dimensions (W x H x D)	76 x 74 x 44 (± 1 mm)			
Cable length	1,800mm (± 10 mm)			
Adapters available	USA	Europe	U.K.	Australie, NZ
				
AC input voltage	100 V – 240 V (± 10%)			
AC input current	Max. 250 mA			
AC input frequency	60 Hz / 50 Hz			
DC output current	12 V 1000 mA			
Output protection	Short circuit			
Operating conditions	0°C to 25°C			
Storage conditions	-25°C to +70°C			
Protection class	II			
Unit casing materials	ABS plastic, in accordance with standard UL 94V-0			

3. Technical data relating to the FlowSense® consumable

Specifications	FlowSense®
Dimensions (W x H x D)	Length: 69.3 mm Flow tube diameter: 11.25 mm FlowSense® inlet and outlet connections comply with EN 5356.
Weight	< 20 g
Shelf life	3 years at room temperature <i>Warning:</i> The shelf life of the FlowSense® may be reduced if inappropriate storage/transport conditions are used for more than 48 hours.

#### 4. Technical data related to electromagnetic compatibility (EMC)

EN

Manufacturer statement and guidance – Electromagnetic interference		
EOLife® is designed to operate in an electromagnetic environment as defined below. The user should ensure that the device is used in such an environment.		
Interference measurements	Compliance	Electromagnetic environment – Guidance
Conducted emissions – Mains terminal disturbance voltage according to CISPR11	Class B: 150 kHz to 30 MHz	Equipment suitable for use in households and buildings used for domestic purposes directly connected to a low voltage power supply.
Radiated emissions – Radiation disturbance according to CISPR11	Class B: 30 MHz to 1 GHz	
Harmonic current emissions according to IEC 61000-3-2	Class A 0 kHz – 2 kHz	/
Voltage fluctuations/flicker according to IEC 61000-3-3	0 kHz to 2 kHz	/
Manufacturer statement and guidance – Electromagnetic immunity		
EOLife® is designed to operate in an electromagnetic environment as defined below. The user should ensure that the device is used in such an environment.		
Immunity tests	IEC 60601 test level	Level of compliance
Electrostatic discharge (ESD) according to IEC 61000-4-2	Contact: ±8kV Air: ±2kV ±4kV, ±8kV, ±15kV	Contact: ±8 kV Air: ±15 kV
Electrical fast transients/bursts according to IEC 61000-4-4	Access by power AC/DC: ± 2 kV, f = 100 kHz	5/50 ns, 100 kHz; ± 2 kV
Shock waves/surges according to IEC 61000-4-5	Access by AC/DC power ports: ± 0.5 kV, ±1 kV (Line/Line)	1.2/50 (8/20) µs LtL: ±1.0 kV
Voltage dips, short interruptions and voltage variations according to IEC 61000-4-11	Min. & max. voltage: 0% UT for 0.5 cycle 0% UT for 1 cycle 70% UT for 25/30 cycles 0% UT for 250/300 cycles	0% UT for 0.5 cycle (1 phase) 0% UT for 1 cycle 70% UT for 25/30 cycles 0% UT for 250/300 cycles (50/60 Hz)
Conducted disturbances according to IEC 61000-4-6	Access by AC/DC power ports: 3 V (0.15 MHz to 80 MHz) 6 V (+ ISM bands) 80% MA, 1 kHz	150 kHz – 80 MHz 3 V ISM/amateur bands 6 V 80% / 1 kHz other as defined in RMF
Power frequency magnetic field according to IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 45 50 Hz or 60 Hz
Radiated disturbances according to IEC 61000-4-3	80 MHz to 2.7 GHz AM modulation 80%, f = 1 kHz 10 V/m: Home	80 MHz to 2.7 GHz Home Healthcare 10 V/m Prof. Healthcare 3 V/m 80% / 1 kHz other as defined in RMF
EOLife® has been tested for immunity to the radio services listed below. If the field strength measured on the site where EOLife® is used exceeds the level of compliance below, it is advisable to monitor EOLife® to ensure that it is functioning properly. If this is not the case, it may be necessary to take additional measures, such as changing the direction or location of EOLife®. Specific examples of RF transmitters known to cause electromagnetic interference are diathermy and electrocautery devices, RFID devices and security systems (such as electromagnetic anti-theft devices and metal detectors). Some of these RF transmitters (such as RFID devices) may be concealed in the area where the device is to be used, potentially exposing it to these RF transmitter fields.		

Proximity fields from RF wireless communications equipment according to IEC 61000-4-3	380 to 390 MHz 27 V/m; PM 50%; 18 Hz	380 to 390 MHz 27 V/m; PM 50%; 18 Hz 430 to 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz 704 to 787 MHz 9 V/m; PM 50%; 217 Hz 800 to 960 MHz 28 V/m; PM 50%; 18 Hz 1700 to 1990 MHz 28 V/m; PM 50%; 217 Hz 2400 to 2570 MHz 28 V/m; PM 50%; 217 Hz 5100 to 5800 MHz 9 V/m; PM 50%; 217 Hz
	430 to 470 MHz 28 V/m; (FM ± 5 kHz, 1 kHz sin) PM; 18 Hz	
	704 to 787 MHz 9 V/m; PM 50%; 217 Hz	
	800 to 960 MHz 28 V/m; PM 50%; 18 Hz	
	1700 to 1990 MHz 28 V/m; PM 50%; 217 Hz	
	2400 to 2570 MHz 28 V/m; PM 50%; 217 Hz	
5100 to 5800 MHz 9 V/m; PM 50%; 217 Hz		

## X. What to do in the event of an incident

If there is a persistent problem or a major adverse event occurs when using EOLife® (e.g. battery problem or sudden device shutdown), ask your approved distributor what procedure to follow and/or contact the manufacturer ARCHEON at this address: [product-request@archeon-medical.com](mailto:product-request@archeon-medical.com).

## XI. Warranty and limits of liability

1. The manufacturer guarantees that the EOLife® device was manufactured in accordance with the technical specifications, good manufacturing practices and other industrial standards and applicable regulations.
2. The manufacturer undertakes to replace or refund any EOLife® product with hidden defects before the warranty expires, providing that it receives the lot/serial number of the defective product.
3. This warranty supersedes any other written or oral, express or implied, statutory or other guarantee, and no commercial or other warranty that differs from this one will be applicable. The only legal remedy in the event of manufacturing defects is the warranty provided above. This warranty does apply in the event of loss, damage, injury or expenses incurred directly or indirectly linked to the use of EOLife®.
4. The manufacturer is released from any responsibility in the event of misuse or improper handling, non-compliance with warnings and instructions, damage occurring after EOLife® is offered for sale, or for any other guarantee given by approved distributors.
5. The manufacturer's warranty is valid for the following periods:

Product	Warranty period
EOLife® electronic unit	2 years
Battery	2 years
Charger	2 years

6. The manufacturer is ARCHEON, 2 Chemin des Aiguillettes, 25000 Besançon, France.

## XII. CE declaration of conformity for medical devices

ARCHEON hereby declares that the product complies with the relevant provisions of Regulation (EU) 2017/745 on medical devices.

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