# EOLifeX

INSTRUCTIONS FOR USE  
INSTRUCTION D'UTILISATION

ARC  EON  
MEDICAL

Patient's type

55%

Child  
(<140cm)

Adult  
(>140cm)

EOlifeX

# EOLifeX

---

[EN] Instructions for use . . . . .	4
[FR] Instruction d'utilisation . . . . .	24

# TABLE OF CONTENTS

EN

- I. Symbols used on product labelling . . . . . 5
- II. Acronyms. . . . . 5
- III. Conditions of use . . . . . 6
  - 1. Indications for use . . . . . 6
  - 2. Patients. . . . . 6
  - 3. Environment of use . . . . . 6
  - 4. Users. . . . . 6
  - 5. Contraindications . . . . . 6
- IV. Safety. . . . . 6
- V. Product description . . . . . 7
- VI. Equipment and accessories. . . . . 8
- VII. Quick start guide of the device . . . . . 9
- VIII. Detailed product description and operating instructions . . . . . 11
  - 1. Using the device. . . . . 11
    - a) Assembly/disassembly . . . . . 11
    - b) EOlife X blind ventilation training mode . . . . . 12
    - c) EOlife X monitored ventilation training mode . . . . . 12
    - d) Visual signals conditions . . . . . 15
    - e) Visual signals prioritisation . . . . . 17
    - f) Automatic shutdown . . . . . 18
  - 2. EOlife X device set-up screen . . . . . 18
    - a) Set-up screen display . . . . . 18
    - b) Device settings. . . . . 19
    - c) Data transfer . . . . . 19
    - d) Update . . . . . 20
    - e) Information . . . . . 21
  - 3. Maintenance. . . . . 21
    - a) Charging the battery . . . . . 21
    - b) Disposal . . . . . 22
- IX. Technical data . . . . . 23
  - 1. Technical data related to the device . . . . . 23
  - 2. Technical data relating to the power supply . . . . . 24
  - 3. Technical data relating to the FlowSense X . . . . . 25
- X. What to do in the event of an incident . . . . . 25
- XI. Warranty and limits of liability . . . . . 25
- XII. Certifications . . . . . 26

The EOlife X equipment is manufactured by ARCHEON,  
2 Chemin des Aiguillettes, 25000 Besançon, France.

## I. Symbols used on product labelling

EN

Symbol	Definition
	Manufacturer's lot number
	Part number
	Serial number
	Quantity
	CE marking
	Warning
	Manufacturer
	Use by date
<b>IP44</b>	Ingress protection against : <ul style="list-style-type: none"><li>• Intrusion of solid objects larger than 1 mm</li><li>• Water splashing from any direction</li></ul>
	Temperature limitation
	Humidity limitation
	Do not dispose of with general waste
<b>UK CA</b>	UKCA UK Conformity Assessed

## II. Acronyms

<b>AHA</b>	American Heart Association
<b>CPR</b>	Cardiopulmonary resuscitation
<b>ERC</b>	European Resuscitation Council
<b>LED</b>	Light-Emitting Diode
<b>WEEE</b>	Waste Electrical and Electronic Equipment

## III. Conditions of use

---

EN

### 1. Indications for use

EOLife X is a radio equipment intended to the training of users on a mannikin to the manual ventilation during cardiopulmonary resuscitation (CPR).

### 2. Patients

EOLife X is indicated to be used on:

- adult simulated patient between 4'7" (140 cm) and 6'7" (200 cm)
- pediatric simulated patient between 1 year (above 2'5" = 75 cm) and 12 years (below 4'7" = 140 cm)

### 3. Environment of use

EOLife X is intended to be used in training session.

### 4. Users

EOLife X is intended to be used by healthcare professionals and rescuers.

### 5. Contraindications

EOLife X is not approved for the following uses:

- In a hyperbaric chamber
- Near magnetic resonance imaging machines

## IV. Safety

---

### Precautions

- o Before using the product, read these instructions carefully. Failure to follow the instructions can lead to improper use and expose the user or persons nearby to risks.
  - Use the device as intended (see "III. Conditions of use").
  - Comply with the contraindications (see "III Conditions of use").
  - Carefully read the safety instructions in these instructions for use.
  - Comply with all the sections in these instructions.
- o The use of defective devices may put the user at risk.
  - Only use the EOLife X device and its accessories if they do not show any visible damage.
- o The EOLife X device, its battery, its charger and the FlowSense X sensor must be stored and used as described in the instructions for use (see "X. Technical data").
- o The device can be safely used outdoors, as it is rainproof, splashproof and cannot be penetrated by solid objects larger than 1 mm (IP44 rating). However, it should be handled with care and kept away from powerful water jets and dust.
- o Always check the battery level before using EOLife X. We advise recharging the battery after each use. Change the battery if the charge is below 20%. We also advise having a spare operational battery in case this is needed.
- o Using a battery with a low charge at temperatures below 0°C greatly reduces the run time and can cause EOLife X to stop functioning after a short time. If the device is used at low temperatures, the battery should be fully charged.
- o Do not remove the battery when EOLife X is switched on.


## Warnings

- ⚠ **EOLife X is not intended for use on humans. EOLife X should only be used on a manikin or patient simulator.**
- ⚠ Avoid using EOLife X next to or in combination with other devices (battery, charger or other) not recommended by ARCHEON as this could cause it to malfunction. If devices must be used next to each other or together, monitor both devices to check that they are functioning properly.
- ⚠ To ensure that EOLife X functions effectively you must follow the instructions provided by ARCHEON.
- ⚠ Never use the EOLife X device or its accessories in a hyperbaric chamber as this could cause an explosion.
- ⚠ The use of accessories, transducers and cables other than those specified or supplied could increase the electromagnetic emissions or decrease the electromagnetic immunity of the device and cause it to malfunction.
- ⚠ Portable high frequency communication devices (e.g. wireless devices, antennas and their cables) in close proximity to EOLife X can interfere with the device and reduce its performance. Keep EOLife X at least 30 cm from portable high frequency communication devices.
- ⚠ Handling the battery improperly could seriously injure the user and persons nearby. The battery must not be disposed of in a fire, opened, deformed or short-circuited. The battery must be protected from humidity, high temperatures and high pressure (see "X. Technical data").
- ⚠ The EOLife X electronic unit, the FlowSense X sensor, the battery and the charger must be discarded at the end of their lifetime as instructed in "IX. 3 d) Disposal".

## V. Product description

- The device is a portable machine that includes the EOLife X electronic unit, its removable and rechargeable battery, its charger and the FlowSense X sensor.


- EOLife X is used with the FlowSense X sensor which is positioned between a non-invasive ventilation interface (mask), invasive ventilation interface (endotracheal tube), or any type of supraglottic airway device, and a standard self-inflating bag to manually ventilate a manikin.

#### Note

The ventilation interface and the self-inflating bag are accessories not supplied by ARCHEON.

- EOLife X allows the user to choose the patient type (adult or child), the ventilation mode (continuous ventilation, CPR with a 30:2 ventilation ratio or CPR with a 15:2 ventilation ratio), the patient's height or the patient's age and the training mode (blind ventilation or monitored ventilation).
  - **Monitored ventilation training:**  
Based on the inspiratory and expiratory flow rates measured during ventilation, EOLife X calculates and records the main ventilation parameters (insufflation volumes, tidal volumes, ventilation frequency and leakage) and gives real-time feedback on the ventilation performed. EOLife X also shows the target ventilation values (in accordance with the guidelines issued by the AHA and ERC international scientific committees) to guarantee adequate ventilation.
  - **Blind ventilation training:**  
Based on the inspiratory and expiratory flow rates measured during ventilation, EOLife X calculates and records the main ventilation parameters (insufflation volumes, tidal volumes, ventilation frequency and leakage) but does not give the user real-time feedback on the ventilation performed and does not show target ventilation values.
- EOLife X records the ventilation data, which can be exported to another device for research and analysis. EOLife X includes a Bluetooth low energy wireless module to communicate with such devices.

## VI. Equipment and accessories

Equipment	Part number
EOLife X	A0000089
EOLife X battery	A0000110
Charger and 4 adapters	A0000029
FlowSense X (x1)	A0000093
FlowSense X (x10)	A0000095

The equipment and accessories can be ordered from approved distributors.


## VII. Quick start guide of the device

This section gives a short step-by-step description of how to use the EOLife X device. There is a more detailed description in the next section, "VIII. Detailed product description and operating instructions".


Assemble EOLife X, attaching the FlowSense® single use sensor and other devices to ventilate the patient (ventilation interface and bag)

EN


- 1** EOLife X assembly with the FlowSense® sensor  
See "VIII 1 a) Assembly"

- 2** Assembly with the ventilation interface

- 3** Assembly with the bag

- 4** Start up EOLife X by pressing the physical ON/OFF button

- 5** Press on the type of the patient :


- ▶ **Child** : the patient is a child under 140 cm (4'7") tall.
- ▶ **Adult** : the patient is an adult above 140 cm (4'7") tall.

- 6** According the patient type (child or adult), press on :


▶ **Age of the child patient**

- purple : 1-2 years
- Yellow : 2-3 years
- White : 3-4 years
- Blue : 5-6 years
- Orange : 7-9 years
- Green : 10-11 years


▶ **The height of the adult patient :**


- small : the patient's height is between 140 cm (4'7") and 160 cm (5'2")
- medium: the patient's height is between 160 cm (5'3") and 180 cm (5'10")
- tall : the patient's height is between 180 cm (5'11") and 200 cm (6'7")


- 7** According to the patient type (child or adult), choose the type of ventilation performed :
- Child :**
- 15:2 mode : alternating between 15 chest compressions and 2 ventilations
  - Continuous mode : continuous ventilations at a frequency recommended by AHA and ERC guidelines


- Adult :**
- 30:2 mode : alternating between 30 chest compressions and 2 ventilation
  - Continuous mode : ventilation at a frequency of 10 cycles per minute

- 8** Choose the type of training :
- Blind ventilation
  - Monitored ventilation

- 9** Ventilation screen
- See VIII 1 b) EOLife X ventilation method
  - See VIII 1 d) Visual signals

Information area (advanced interface)


- Vi = insufflated volume
- Vt = tidal volume
- Freq. = ventilation frequency

Bar graph (real-time Vi)


Visual signal display area

Main screen (advanced interface)

- 10** Display the summary screen by pressing the physical ON/OFF button


- 11** Select the date, time and language and choose to deactivate the leakage alarm and the followed guideline

- 12 Switch off Eolife X by pressing the ON/OFF button for 3 seconds or by quick pressing the ON/OFF icon on the touch screen

Switching off using the summary screen

- 13 Disconnect the ventilation bag and interface (mask or endotracheal tube) and unclip the FlowSense X sensor by pulling the tab outwards, as shown by the arrow

Disconnecting Eolife X

## VIII. Detailed product description and operating instructions

Eolife X is controlled via the touch screen and the physical ON/OFF button.

 **Caution**

Avoid using any object that could scratch the screen.

### 1. Using the device

#### a) Assembly/disassembly


##### Assembly

FlowSense X is factory calibrated. Once removed from its packaging FlowSense X is ready to be used. Remove the FlowSense X sensor from its packaging by using the notch to tear the antistatic bag to avoid damaging the sensor.

Connect the FlowSense X sensor to the Eolife X unit, taking care not to expose the printed circuit board to dust or splash it with water.


First insert the part of the FlowSense X with a notch into the connector ridge, and then press firmly on the sensor so that it comes into contact with the connector pins. You will hear a click when the connector tab holds the FlowSense X sensor in place.

Check that the FlowSense X sensor is held properly in place before starting up the Eolife X unit.

Connect the ventilation interface (mask or endotracheal tube) to the compatible end of the FlowSense X sensor. The sensor ends are shaped differently to ensure that the interface cannot be connected to the wrong end.

Connect the bag to the other end of the FlowSense X sensor. The device must be assembled as shown here:


## Disassembly

To remove the accessories from the product, take the FlowSense X sensor out of the ventilation interface (mask or endotracheal tube) by gently pulling vertically on the FlowSense X sensor. Then remove the ventilation bag at the other end in the same way.

### b) EOLife X blind ventilation training mode

In blind ventilation mode, EOLife X calculates the volumes insufflated and expired, the ventilation frequency and leakage, estimates the tidal volumes, and records the data. After recording the ventilation session, EOLife X provides a statistical report on the quality of the ventilation performed.

### c) EOLife X monitored ventilation training mode

In monitored ventilation mode, EOLife X gives the user real-time feedback on the quality of the ventilation performed based on a cyclic analysis of inspiratory and expiratory flow rates. It must therefore be positioned proximally to be able to measure both the inspiratory and expiratory flow rates, calculate the volumes insufflated and expired, the ventilation frequency and leakage, estimate the tidal volumes, and generate user information signals and visual signals if the ventilation parameters are unacceptable.


EOLife X also uses algorithms to filter out small airflows generated by chest compressions so as not to display anomalous volume and ventilation frequency values and to avoid measurement artefacts. During ventilation, be sure to generate enough flow to be detected by the device.

#### **Warning**

When ventilating with a mask, ensure it is correctly held on the manikin's face during insufflation and for the entire duration of expiration so that EOLife X can measure the volume expired.

**Note**  
 Check that the manikin's airways are completely sealed so that all of the air insufflated via the mouth is fully expired. If this is not the case, you have to deactivate the leakage alarm in the EOlife X settings (see VIII 2 b)). EOlife X will be able to analyse the ventilation cycles or calculate the ventilation parameters correctly by simulating the expiration.

The main screen guides the user by displaying the ventilation parameters and visual guide showing the ideal air/oxygen volume to administer and when to administer it. The main screen also uses visual signals to warn the user if the manikin is not ventilated correctly.

1. Displays the patient's height or the patient's age selected and quick access key to change the patient's height or the patient's age
2. Displays the chosen ventilation mode and quick access key to change the ventilation mode.
3. Rotation icon and key to rotate the screen 180° to adjust the direction of the screen to the user's position during ventilation.  
 The user can perform ventilation and change the patient's height, the patient's age or the ventilation mode. The screen returns to its normal position if the user:
  - Presses again on the rotation icon
  - Presses the physical ON/OFF button to show the summary screen
  - Presses the physical ON/OFF button for 3 seconds to switch the device off
4. Displays the battery level (between 0% and 100%).
  - o See "VIII. 3 a) Charging the battery".
5. Displays the actual volume insufflated (Vi) for each ventilation cycle (value in mL).
6. Displays the actual tidal volume (Vt) (for each of ventilation cycle (value in mL).


**Note**  
 The tidal volume is the estimated volume of air/oxygen actually reaching the patient's lungs based on both the measurement of the volume of air expired and the leak calculated during the insufflation/expiration phases. The ERC and AHA recommend that adult and child patient in cardiopulmonary arrest should be ventilated with a tidal volume of 6 to 8 mL/kg of Ideal body weight.

7. Displays the trend value of the ventilation frequency (Freq.) (value in cycles/minute) taking into account the frequency of the last ventilation cycles performed :  
 For the 30:2 and the 15:2 mode : value displayed from the third ventilation cycle.  
 For the continuous mode : value displayed from the fourth ventilation cycle.


**Note**

Adult (ERC and AHA)	10 ventilations per min
Child from 1 to 6y (ERC)	20 ventilations per min
Child from 7 to 11y (ERC)	15 ventilations per min
Child from 1 to 11y (AHA)	25 ventilations per min

- 8. Real-time indication of how to ventilate:
  - ▶ A bar graph shows whether the volume insufflated/exhaled is adequate:

- When the user ventilates, the bar graph shows in real time the quantity of air/oxygen that must be insufflated based on ERC and AHA guidelines.
  - ▶ Display in continuous mode:

- In continuous resuscitation mode, EOLife X generates a visual countdown to show the user when to ventilate to deliver the 'target frequency' breaths per minute. When the green lung is displayed, the user must ventilate.
  - Display in 30:2 and 15:2 mode:

In 30:2 mode, the countdown function is disabled as the person ventilating must wait for 30 chest compressions before delivering 2 ventilations.  
In 15:2 mode, the countdown function is disabled as the person ventilating must wait for 15 chest compressions before delivering 2 ventilations.  
In this mode the ventilation frequency depends on the person (or machine) performing the chest compressions.  
The green lung is shown permanently to indicate that the device is waiting to detect a ventilation to display the bar graph.

### 9. Visual signal display area :

- Good Ventilation

all the ventilation parameters calculated are within the appropriate range of values as defined in the international guidelines.
- Insufficient Volume

the volume of air insufflated to the patient is significantly lower than the values given in the international guidelines.
- Excessive Volume

the volume of air insufflated to the patient is significantly higher than the values given in the international guidelines.
- Low Frequency

the ventilation frequency is much significantly than the values given in the international guidelines.
- High Frequency

the ventilation frequency is significantly higher than the values given in the international guidelines.
- Important Leakage

the difference calculated between the volume of air/oxygen insufflated and the volume exhaled is too great. This may be because the mask is not positioned correctly, due to improper intubation (oesophageal intubation) or because the ventilation circuit is not properly sealed.

#### Note

See details of visual signals in "VIII.1 c) Visual signals".

### d) Visual signals conditions

#### Visual messages related to ventilation parameters

Visual messages may appear during ventilation when the ventilation parameters exceed the threshold values and are persistently outside the ranges given below.

The screen shows the visual messages in the visual signal display area (lower part of the screen) as a message.

#### The visual messages concern:

- **The tidal volume (Vt):**

Visual signal conditions	Cause		Solution
	Patient height	Tidal volume value	
Insufficient volume	Small	< 225 mL	Insufflate a volume of air within the appropriate range as defined in the international guidelines
	Medium	< 285 mL	
	Tall	< 350 mL	
	1-2 years	< 45 mL	
	2-3 years	< 55 mL	
	3-4 years	< 70 mL	
	5-6 years	< 90 mL	
	7-9 years	< 115 mL	
10-11 years	< 140 mL		

Visual signal conditions	Cause		Solution
	Patient height	Tidal volume value	
Excessive volume	Small	> 430 mL	Insufflate a volume of air within the appropriate range as defined in the international guidelines
	Medium	> 565 mL	
	Tall	> 710 mL	
	1-2 years	> 100 mL	
	2-3 years	> 120 mL	
	3-4 years	> 150 mL	
	5-6 years	> 190 mL	
	7-9 years	> 240 mL	
	10-11 years	> 300 mL	

**Note**

- In continuous ventilation mode the visual signal condition is calculated based on the trend value of the volumes. If the trend value and the current value of the tidal volume are as described in the table above, the tidal volume visual signal is triggered.
- In 30:2 and 15:2 CPR mode the visual signal condition is calculated based on the mean value of the tidal volumes during two consecutive ventilation cycles.

• **The ventilation frequency (Freq.):**

Visual signal conditions	Cause	Ventilation frequency value	Solution
Low frequency (continuous mode)	Adult	No insufflation during 9s	Ventilate at an appropriate frequency, as recommended in the international guidelines
	Child: 1-4y (ERC) and AHA	No insufflation during 6s	
	Child: 5-11y (ERC)	No insufflation during 6.5s	
Low frequency (30:2 mode)	Adult	No inspiration for 24 seconds	
Low frequency (15:2 mode)	Child: 1-4y (ERC) and AHA	No insufflation during 14s	
	Child: 5-11y (ERC)	No insufflation during 14.5s	
High frequency	Adult	> 15 bpm	
	Child: 1-4y (ERC) and AHA	> 30 bpm	
	Child: 5-11y (ERC)	> 25 bpm	

**Note**

- In continuous ventilation mode the frequency is calculated based on the trend value of the ventilation frequency. If the trend value and the current value of the frequency are as described in the table above, the ventilation frequency visual signal is triggered.
- In 30:2 and 15:2 CPR mode the visual signal condition is calculated based on the actual value of the pause time.

• Leakage:

Visual signal conditions	Leakage value	Solution
Important leakage for adult patient	> 40%	Depending on the cause of the leakage, several solutions are possible (non-exhaustive list): <ul style="list-style-type: none"> <li>○ Reposition the tracheal tube if the manikin is improperly intubated (selective intubation, tube not deep enough inside the trachea) or inflate the bag more if the ventilation circuit is not properly sealed</li> <li>○ Reposition the mask properly on the manikin's face to ensure that the ventilation circuit is properly sealed</li> <li>○ Completely free the manikin's airways</li> </ul>
Important leakage for child patient	> 50%	

**Note**

- In continuous ventilation mode leakage is calculated based on the trend value of the leakages. If the trend value and the current leakage value are as described in the table above, the leakage visual signal is triggered.
- In 30:2 and 15:2 CPR mode the visual signal condition is calculated based on the mean value of the leakages during two consecutive ventilation cycles.

**e) Visual signals prioritisation**

All visual signal messages appear in the same display area, therefore a system is necessary to prioritize the visual signals. Visual signals prioritization also enables the user to focus first on correcting the most critical problem in terms of patient safety, and avoids the additional stress of multiple simultaneous visual signals.

The visual signals were prioritized with clinical experts based on the risk level:

**Visual signal message linked to the FlowSense X® sensor**

During ventilation, if the FlowSense X sensor is defective or not properly connected or if a foreign object on the connector interferes with the signal transmission, an alarm screen is displayed:

FlowSense X failure screen

You must reconnect a functional FlowSense X sensor as described above to clear the "FlowSense X failure" screen.

**Note**

The "FlowSense X failure" alarm screen can only be displayed from the main screen.

### EOLife X device failure

If there is a device failure, EOLife X will be unusable: the screen switches off (or no longer comes on) and the red LED in the top right corner of the screen comes on. See "XI. Troubleshooting".

### f) Automatic shutdown


If there is no user action (a press of the physical ON/OFF button, a touch of the screen when the device is running or the detection of the inspiration phase when a ventilation sequence is running,) for 1 hour, the device turns itself off completely.

To turn the device on again, simply press the physical ON/OFF button.

## 2. EOLife X device set-up screen


### a) Set-up screen display

1. Start up EOLife X by pressing the physical ON/OFF button.
2. Press normally on the physical ON/OFF button to show the summary screen.


Summary screen

3. Next, select the [wrench ] to display the EOLife X set-up screen.

Summary screen


Set-up screen

## b) Device settings

### o Leakage alarm settings

When the manikin's airways are not completely sealed and the air insufflated via the mouth is not fully expired, you have to deactivate the leakage alarm in the device settings because EOLife X will not correctly work. By deactivating the leakage alarm, EOLife X will simulate an expiration and will be able to analyse the ventilation cycles or calculate the ventilation parameters correctly.


EN

### o Frequency settings

For a child patient, AHA and ERC guidelines do not have the same recommendations about the ventilation frequency. EOLife X offers the possibility to choose between following AHA or ERC guidelines.

- AHA guideline for children between 1 year and 11 years : 25 cycles per minutes
- ERC guideline for children between 1 year and 6 years : 20 cycles per minutes
- ERC guideline for children between 7 years and 11 years : 15 cycles per minutes

## c) Data transfer

EOLife X records the ventilation parameters in each ventilation cycle. It can record for an estimated period of 10 hours. This data can then be transferred via Bluetooth to a compatible device using the EOLife Connect® app to extract and analyse the data. To avoid losing the data recorded, it is advisable to transfer the data recorded by EOLife X after each use.

### Note


Transferring data requires the use of a device compatible with EOLife Connect®, not included with the equipment sold with EOLife X. This device can be either a tablet or a smartphone running on one of these operating systems: Android or IOS.

The EOLife Connect® app must be downloaded from the App Store or Google Play Store and installed on your device to be able to use the EOLife X data transfer feature. Please contact an approved distributor or the manufacturer for more information.

### Note

When data is being transferred from the device, the EOLife X screen is frozen and the red LED at the top left of the screen may come on. Once the transfer is complete, the LED will turn off and the device will be operational again.

Download the data recorded by EOLife® from the set-up screen by selecting [the  icon]. The EOLife® Bluetooth connection is automatically turned on and a message appears.

Bluetooth activation screen

Download the data recorded by EOLife X from EOLife Connect®:

1. Start up EOLife Connect®.
2. Press on [the  icon] on the home screen. The EOLife Connect® Bluetooth is enabled and the device searches for nearby EOLife X devices.
3. On the connection screen, select the serial number of the EOLife X device you want to connect with to perform the data transfer.
4. A circular progress bar shows the progress of the data transfer. A confirmation message appears once the data transfer is complete.
5. Press on one of the two options to view either the last downloaded ventilation or all downloaded ventilations.

## d) Update

EOLife X has software that can be updated via Bluetooth using the EOLife Connect® app. Your approved distributor or the manufacturer will notify you when an update is available.

### Note

This update requires the use of a device compatible with EOLife Connect®, not included with the equipment sold with EOLife X. This device can be either a tablet or a smartphone running on one of the following operating systems: Android or IOS.


The EOLife Connect® app must be downloaded from the App Store or Google Play Store and installed on your device to be able to use the EOLife X software update feature.

Please contact an approved distributor or the manufacturer for more information.

### Note

During the device update, the EOLife X screen is off and the red LED at the top left of the screen is on. Once the update is complete, the LED will turn off and the device will start up again.

Turn on the EOLife X's Bluetooth from the set-up screen by selecting [the  icon]. It turns on automatically and a message appears.

Bluetooth activation screen

Update the EOLife X software from EOLife Connect®:

1. Start up EOLife Connect®.
2. Press on [the  icon] on the home screen. The EOLife Connect® Bluetooth is enabled and the device searches for nearby EOLife X devices.
3. On the update screen, select the serial number of the EOLife X device you want to connect with to perform the update.
4. A progress bar shows the progress of the update. A confirmation message appears once the update is complete.

## e) Information

The information contained in the EOLife X device includes:

- o The EOLife X serial number
- o The software version
- o The operating conditions
- o Temporary operating conditions
- o The list of certifications
- o Storage conditions
- o The manufacturer
- o The manufacturer's address

## 3. Maintenance

### a) Charging the battery


Recharging the EOLife X battery:

**Note**

The battery must be removed from the EOLife X unit before recharging.  
Turn off the device before battery removal.

1. Remove the battery.
2. Plug the jack on the charger cable provided into the battery.

3. Plug the charger cable to the mains socket using the appropriate adapter for your country.

Adapter name	Male adapter B	Male adapter E	Male adapter G	Male adapter I
Image				
Compatible sockets	Type A Type B	Type C Type E Type F	Type G	Type I

**Warnings**

- o Only recharge the EOLife X battery using the charger provided.
- o Only use the EOLife X battery referred to on the list of accessories and equipment (part number not available).
- o If the charge is too low, do not interrupt ventilation to replace the battery. Continue ventilation without real-time feedback from EOLife X.

Always check the battery level before placing the device in the emergency bag. The battery level is shown in the top right corner of the start-up screen.

The battery icon fills as the charge increases and also changes colour:

	Green Run time ≥ 50%	Orange Run time 20%-49%	Red Run time < 20%
Charge			

 **Caution**

When the  symbol appears on the screen, the battery is probably defective and you should use another battery to avoid a power failure during ventilation.

**b) Disposal**

At the end of the EOlife X device's lifetime (lifetime of the various parts specified in "X. Technical data"), it should be disposed of as follows:

The EOlife X electronic unit, the FlowSense X sensor, battery and charger contain electronic components and must be disposed of in an appropriate recycling facility, in accordance with European Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). They must not be disposed of with household waste. Take them to an electrical and electronic equipment recycling collection point. They must be disposed of in accordance with local environmental regulations on waste disposal. For more detailed information on the handling, collection and recycling of this product, contact your environmental officer or local council.

## IX. Technical data

### 1. Technical data related to the device




Specifications	Eolife X
Dimensions (W x H x D)	130 mm x 75 mm x 30 mm (5.11 inches x 2.95 inches x 1.18 inches)
Weight	170g ± 5 g (5.3 oz ± 0.011 oz)
Operating conditions	<ul style="list-style-type: none"><li>• Temperature between 0°C and 40°C (32°F and 104°F)</li><li>• Relative humidity between 15% and 95% (non-condensing)</li><li>• Atmospheric pressure range of 620 hPa (altitude of 4000 m) to 1,060 hPa (altitude of -500 m)</li></ul>
Storage/transport conditions	<ul style="list-style-type: none"><li>• Temperature between -10°C and 35°C (-14°F and 95°F)</li><li>• Relative humidity between 10% and 90% (non-condensing)</li></ul>
Lifetime	5 years
Run time	5 hours
IP rating for protection against solids, dust and ingress of water	IP44 (in use configuration, i.e., Eolife X, its battery and FlowSense X connected)
Electromagnetic compatibility (EMC) according to EN 60601-1-2	The control parameters and threshold values can be obtained from the manufacturer.
BLE module frequency range	2,4 GHz.
Signal strength	Max. 4 dBm.
Certification	Directive 2014/53/EU, FCC, IC, MIC, KC.
Screen	2.4 inches. Resolution 320 x 240 pixels.
Measurement accuracy	<ul style="list-style-type: none"><li>• Vi (volume insufflated): ± 8.8% of the actual value measured under normal conditions of use</li><li>• Vt (tidal volume) without leakage: ± 8.8% of the actual value measured under normal conditions of use</li><li>• Freq (ventilation frequency): ± 1 cycle per minute</li></ul>

EN

## 2. Technical data relating to the power supply

EN

Specifications	EOLife X battery
Dimensions (W x H x D)	73 mm x 40 mm x 14 mm
Weight	50 g ± 5 g
Shelf life	12 months (with a storage temperature between -10°C and 35°C (-14°F and 95°F) Charging the device every 6 months allows to increase battery shelf-life.
Type	Li-ion
Nominal capacity	1,900 mAh (7.03 Wh)
Nominal voltage	3.7 V
Charging voltage	12V.
Charge time (0% to 95%)	Approx. 2 hours
Recommended charging temperature	+10°C to 45°C
Certifications	IEC 62133 UN38.3 UL1642

Specifications	Charger								
Weight	206 g including accessories								
Unit dimensions (L x H x D)	76 x 74 x 44 (± 1 mm)								
Cable length	1,800 mm (± 10 mm)								
Adapters available	<table border="1"> <thead> <tr> <th>USA</th> <th>Europe</th> <th>UK</th> <th>Australia, NZ</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	USA	Europe	UK	Australia, NZ				
USA	Europe	UK	Australia, NZ						
									
AC input voltage	100 V – 240 V (± 10%)								
AC input current	Max. 250 mA								
AC input frequency	60 Hz / 50 Hz								
DC output current	12 V - 1A								
Output protection	Short circuit								
Operating conditions	0°C to +25°C								
Storage conditions	-25°C to +70°C								
Protection class	II								
Unit casing materials	ABS plastic, in accordance with standard UL 94V-0								

### 3. Technical data relating to the FlowSense X

Specifications	FlowSense X
Dimensions (W x H x D)	Length: 69.3 mm Flow tube diameter: 11.25 mm FlowSense X inlet and outlet connections comply with EN 5356.
Weight	< 20 g
Shelf life	3 years at room temperature

EN

## X. What to do in the event of an incident

If there is a persistent problem or a major adverse event occurs when using EOLife X (e.g. battery problem or sudden device failure), ask your approved distributor what procedure to follow and/or contact the manufacturer ARCHEON at this address: [product-request@archeon-medical.com](mailto:product-request@archeon-medical.com).

## XI. Warranty and limits of liability

1. The manufacturer guarantees that the EOLife X device was manufactured in accordance with the technical specifications, good manufacturing practices and other industrial standards and applicable regulations.
2. The manufacturer undertakes to replace or refund any EOLife X product with hidden defects before the warranty expires, providing that it receives the lot/serial number of the defective product.
3. This warranty supersedes any other written or oral, express or implied, statutory or other guarantee, and no commercial or other warranty that differs from this one will be applicable. The only legal remedy in the event of manufacturing defects is the warranty provided above. This warranty does apply in the event of loss, damage, injury or expenses incurred directly or indirectly linked to the use of EOLife X.
4. The manufacturer is released from any responsibility in the event of misuse or improper handling, non-compliance with warnings and instructions, damage occurring after EOLife X is offered for sale, or for any other guarantee given by approved distributors.
5. The manufacturer's warranty is valid for the following periods:

Product	Warranty period
EOLife X electronic unit	2 years
Battery	2 years
Charger	2 years

6. The manufacturer is ARCHEON, 2 Chemin des Aiguillettes, 25000 Besançon, France.

## XII. Certifications

EN

	<p>ARCHEON hereby declares that the product complies with the requirements of Council Directive 2014/53/EU on radio equipment (RED) and Council Directive 2011/65/EU on the restriction of the use of certain hazardous substances. The full text of the declaration of conformity is available on the internet at <a href="http://www.archeon-medical.com">www.archeon-medical.com</a>.</p>
	<p>ARCHEON hereby declares, that the product complies with the requirements of the Radio Equipment Regulations 2017/1206 (RER) applicable in Great Britain. The full text of the declaration of conformity is available on the internet at <a href="http://www.archeon-medical.com">www.archeon-medical.com</a>.</p>
<p><b>Federal Communications Commission</b></p> <p><b>Statement and Industry Canada Statements</b></p>	<p>This device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:</p> <ol style="list-style-type: none"> <li>1. The device may not cause harmful interference.</li> <li>2. This device must accept any interference received, including interference that may cause undesired operation.</li> </ol> <p>Caution: Changes or modifications not expressly approved by ARCHEON could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:</p> <ul style="list-style-type: none"> <li>• Reorient or relocate the receiving antenna.</li> <li>• Increase the separation between the equipment and receiver.</li> <li>• Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.</li> <li>• Consult the dealer or an experienced radio/TV technician for help.</li> </ul> <p>Contains FCC ID: SQGBL652 Contains IC: 3147A-BL652</p>
<p>Certificate of Radio Equipment in Japan (MIC)</p>	<p>No: 201-160415/01</p>   201-160415
<p>Certificate of Broadcasting and Communication Equipments (Korea)</p>	<p>R-C-LAI-BL652-SA</p>

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## XII. Certifications

	<p>ARCHEON déclare par la présente que le produit est conforme aux exigences essentielles de la Directive du Conseil 2014/53/UE relative aux équipements radioélectrique (RED) et de la Directive du Conseil 2011/65/UE relative à la limitation de l'utilisation de certaines substances dangereuses.</p> <p>Le texte intégral de la déclaration de conformité est disponible sur internet à l'adresse suivante <a href="http://www.archeon-medical.com">www.archeon-medical.com</a>.</p>
	<p>ARCHEON déclare par la présente que le produit est conforme aux exigences du Règlement 2017/1206 sur les équipements radio (RER) applicable en Grande Bretagne.</p> <p>Le texte intégral de la déclaration de conformité est disponible sur internet à l'adresse suivante <a href="http://www.archeon-medical.com">www.archeon-medical.com</a>.</p>
<p><b>Federal Communications Commission</b></p> <p><b>Statement and Industry Canada Statements</b></p>	<p>This device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:</p> <ol style="list-style-type: none"> <li>1. The device may not cause harmful interference.</li> <li>2. This device must accept any interference received, including interference that may cause undesired operation.</li> </ol> <p>Caution: Changes or modifications not expressly approved by ARCHEON could void the user's authority to operate the equipment.</p> <p>This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:</p> <ul style="list-style-type: none"> <li>• Reorient or relocate the receiving antenna.</li> <li>• Increase the separation between the equipment and receiver.</li> <li>• Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.</li> <li>• Consult the dealer or an experienced radio/TV technician for help.</li> </ul> <p>Contains FCC ID: SQGBL652 Contains IC: 3147A-BL652</p>
<p><b>Certificate of Radio Equipment in Japan (MIC)</b></p>	<p>No: 201-160415/01</p>   201-160415
<p><b>Certificate of Broadcasting and Communication Equipments (Korea)</b></p>	<p>R-C-LAI-BL652-SA</p>

FR

A0000103-Rev02-22AUG2023